

SMALL CHIP... BIG RETURN

WHY MICROCHIP?

- Sadly 1 in 3 pets will go missing in their lifetime.
- Microchips are the most trusted way of identifying a pet.
- A collar tag is not enough. If your dog's collar breaks or its collar tag falls off or becomes illegible, a microchip provides a permanent, tamper proof means of animal identification.
- Pets can be quickly matched to their current owner, minimizing the expense of housing and feeding while waiting to be returned or re-homed.
- Animal control officers are trained and equipped to scan animals found.
- The implantable microchip provides an electronic means of identification.

We speak for pets everywhere!

Pets with microchips are up to 20 times more likely to be reunited with their owners.



LOST PET... FOUND FAST

**"KNOW YOUR PET IS PROTECTED
24-7-365"**

- It is vital to keep your contact information up to date so that you can be reached.
- Please call us or visit our webpage to update your information.

ELECTRONIC PAWPRINT

Just like a fingerprint no two microchips have the same ID number.

**JUST VISIT US ONLINE @
eidap.com**

or fax your completed form to us
1-866-467-5160

You can call us at
1-888-346-8899

or by email
info@eidap.com

ASK YOUR VETERINARIAN about EIDAP
injectable microchip identification.



**POINTING THE
WAY HOME**

INJECTABLE MICROCHIP IDENTIFICATION

for dogs, cats, equines & all
other species



A BETTER PET TAG

WHY CHOOSE EIDAP TO IDENTIFY YOUR PET?

- Our Trovan® microchips are completely manufactured in Germany to the highest quality standards with over 280 million transponder products sold globally over a period of 25 years!
- When your pet is micro chipped with an EIDAP Inc. microchip, they become a member of the largest and quickest accessible database for animal recovery.
- The EIDAP animal registry is available 24 hours a day and 365 days a year through the internet or toll free, providing all necessary information needed to reunite a lost pet with its owner.
- EIDAP Microchips conform to all international standards (ISO) and can be read worldwide.
- Registration information can be updated at any time – free of charge.
- Microchip numbers are registered in the secure database prior to shipment from EIDAP animal registry to your veterinarian.
- EIDAP's transponders are encapsulated using biocompatible glass, which is impermeable to fluids and non-toxic in the animal.

COMMON QUESTIONS... ANSWERED

- The microchip is placed between your Pet's shoulder blades under a veterinarian's supervision.
- Implantation is quick, easy and virtually painless - similar to a vaccine injection - and can be performed during a regular veterinary check-up visit.
- After implantation, a small layer of connective tissue forms around the microchip, preventing migration.
- The microchip has no power supply, battery or moving parts. It is designed with an operating life of over 20 years and is guaranteed for the life of the animal.

DID YOU KNOW?

- You can upgrade your basic registration to include medical alerts, secondary contact information and even a photo of your pet for added piece of mind!
- It is important to note that this is not a tracking microchip that can be used to pinpoint an animal's location.
- **EIDAP & TROVAN offer multiple RFID solutions to some of the world's largest corporations and most prestigious Universities.**

It's no wonder 300 Zoos worldwide, 150 government agencies and the WWF chose us!



EIDAP & TROVAN® ARE THE OFFICIAL MICROCHIP OF THE CANADIAN KENNEL CLUB (CKC) & THE AMERICAN KENNEL CLUB (AKC)